

Nucleolin Rabbit mAb

ABZ3343


Product Overview

Name	Nucleolin Rabbit mAb
Accession(Primary)	P19338
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 77 kDa; Observed MW: 100 kDa
Recommended Dilution	WB-1:2000-1:5000

Target information (P19338)

Gene Name NCL

About Gene Symbol



Western blot analysis of Nucleolin in HeLa, A549, HL-60 lysates using Nucleolin antibody.

