


Nup107 Rabbit mAb

ABZ3344

Product Overview

Name	Nup107 Rabbit mAb
Accession(Primary)	P57740
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 106 kDa; Observed MW: 106 kDa
Recommended Dilution	WB-1:500-1:1000; IP-1:20-1:50

Target information (P57740)

Gene Name NUP107**About Gene Symbol**

Western blot analysis of Nup107 in A549, HL-60, U251 lysates using NUP107 antibody.

