

NUP155 Rabbit mAb


ABZ3345

Product Overview


Name	NUP155 Rabbit mAb
Accession(Primary)	O75694
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 155 kDa; Observed MW: 155 kDa
Recommended Dilution	WB-1:500-1:1000; ICC/IF-1:500-1:2000; FC-1:100-1:200

Target information (O75694)

Gene Name	NUP155
About Gene Symbol	NUP155 (Nucleoporin 155) is a nuclear pore complex component involved in epigenetic regulation. It is an important research target in epigenetic regulation, chromatin biology, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Immunocytochemistry analysis of NUP155 (green) in HeLa using NUP155 antibody, and DAPI (blue).



Western blot analysis of NUP155 in HeLa, A549, HL-60 lysates using NUP155 antibody.