

NUP50 Rabbit mAb

ABZ3346


Product Overview

Name	NUP50 Rabbit mAb
Accession(Primary)	Q9UKX7
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 50 kDa; Observed MW: 50 kDa
Recommended Dilution	WB-1:1000-1:5000; ICC/IF-1:200-1:500


Target information (Q9UKX7)

Gene Name NUP50

About Gene Symbol



Western blot analysis of NUP50 in HeLa, C6
lysates using NUP50 antibody.



Immunocytochemistry analysis of NUP50 (green)
in hela using NUP50 antibody, and DAPI(blue)