


P2Y13 Rabbit mAb

ABZ3348

Product Overview

Name	P2Y13 Rabbit mAb
Accession(Primary)	Q9BPV8
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 41 kDa; Observed MW: 41 kDa
Recommended Dilution	WB-1:1000-1:5000

Target information (Q9BPV8)

Gene Name P2RY13**About Gene Symbol**

Western blot analysis of P2Y13 in HeLa, A549
lysates using P2Y13 antibody.

