

p57 Kip2 Rabbit mAb


ABZ3349

Product Overview


Name	p57 Kip2 Rabbit mAb
Accession(Primary)	P49918
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 32 kDa; Observed MW: 57 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200; IP-1:10-1:100

Target information (P49918)

Gene Name	CDKN1C
About Gene Symbol	



Western blot analysis of p57 Kip2 in A549, U2OS lysates using p57 Kip2 antibody.



Immunohistochemistry analysis of paraffin-embedded Human breast cancer using p57 Kip2 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.