

PLC gamma 1 Rabbit mAb

ABZ3350



Product Overview

Name	PLC gamma 1 Rabbit mAb
Accession(Primary)	P19174
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 149 kDa; Observed MW: 149 kDa
Recommended Dilution	WB-1:1000-1:5000; ICC/IF-1:50-1:100


Target information (P19174)

Gene Name PLCG1

About Gene Symbol



Immunocytochemistry analysis of Phospho-lipase C gamma 1 (green) in Hela using Phospho-lipase C gamma 1 antibody, and DAPI(blue).



Western blot analysis of Phospho-lipase C gamma 1/PLCgamma1 in U2OS, C6 lysates using Phospho-lipase C gamma 1/PLCgamma1 antibody.