


Pleckstrin Rabbit mAb

ABZ3351

Product Overview

Name	Pleckstrin Rabbit mAb
Accession(Primary)	P08567
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 40 kDa; Observed MW: 40 kDa
Recommended Dilution	WB-1:1000-1:5000

Target information (P08567)

Gene Name PLEK**About Gene Symbol**

Western blot analysis of Pleckstrin in HeLa, HL-60 lysates using Pleckstrin antibody.

