

# Poliovirus Receptor/PVR Rabbit mAb


ABZ3352

## Product Overview


|                             |                                               |
|-----------------------------|-----------------------------------------------|
| <b>Name</b>                 | Poliovirus Receptor/PVR Rabbit mAb            |
| <b>Accession(Primary)</b>   | P15151                                        |
| <b>Host Species</b>         | Rabbit                                        |
| <b>Clonality</b>            | Monoclonal                                    |
| <b>Isotype</b>              | IgG                                           |
| <b>Species Reactivity</b>   | Human                                         |
| <b>Calculated MW</b>        | Calculated MW: 45 kDa; Observed MW: 60-80 kDa |
| <b>Recommended Dilution</b> | WB-1:500-1:1000; IHC-1:200-1:500              |

## Target information (P15151)

|                          |                                                                                                                                                                                                                                          |
|--------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Gene Name</b>         | PVR                                                                                                                                                                                                                                      |
| <b>About Gene Symbol</b> | PVR (PVR Cell Adhesion Molecule) is a membrane receptor involved in adaptive immunity. It is an important research target in infectious disease and is widely studied to elucidate cellular mechanisms and identify therapeutic targets. |



Western blot analysis of Poliovirus Receptor/PVR in HeLa, A549, U2OS lysates using Poliovirus Receptor antibody.



Immunohistochemistry analysis of paraffin-embedded Human Cholangiocarcinoma using Poliovirus Receptor/PVR antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.