

PPCS Rabbit mAb

ABZ3354


Product Overview

Name	PPCS Rabbit mAb
Accession(Primary)	Q9HAB8
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 34 kDa; Observed MW: 34 kDa
Recommended Dilution	WB-1:1000-1:5000

Target information (Q9HAB8)

Gene Name PPCS

About Gene Symbol



Western blot analysis of PPCS in HeLa, Jurkat lysates using PPCS antibody.

