

PRPF4 Rabbit mAb


ABZ3355

Product Overview


Name	PRPF4 Rabbit mAb
Accession(Primary)	O43172
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 58 kDa; Observed MW: 58 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:100-1:200; FC-1:50-1:100

Target information (O43172)


Gene Name	PRPF4
About Gene Symbol	PRPF4 (Pre-mRNA Processing Factor 4) is a protein kinase and signaling regulator involved in translational regulation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of PRPF4 in Raw264.7, C6 lysates using PRPF4 antibody.



Western blot analysis of PRPF4 in HeLa, A549, HL-60 lysates using PRPF4 antibody.



Immunohistochemistry analysis of paraffin-embedded Human Cholangiocarcinoma using PRPF4 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.