

Quiescin Q6 Rabbit mAb


ABZ3356

Product Overview

Name	Quiescin Q6 Rabbit mAb
Accession(Primary)	O00391
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 83 kDa; Observed MW: 83 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:100-1:200; FC-1:50-1:100; IP-1:20-1:50

Target information (O00391)

Gene Name	QSOX1
About Gene Symbol	



Western blot analysis of Quiescin Q6 in A549, U2OS lysates using QSOX1 antibody.