


Rab29 Rabbit mAb

ABZ3357

Product Overview

Name	Rab29 Rabbit mAb
Accession(Primary)	O14966
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 23 kDa; Observed MW: 23 kDa
Recommended Dilution	WB-1:500-1:1000; IP-1:20-1:50

Target information (O14966)

Gene Name RAB29**About Gene Symbol**

Western blot analysis of RAB29 in HeLa, A549, HL-60, U251 lysates using Rab29 antibody.

