

RASA1 Rabbit mAb


ABZ3359

Product Overview


Name	RASA1 Rabbit mAb
Accession(Primary)	P20936
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 116 kDa; Observed MW: 140 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:100

Target information (P20936)

Gene Name	RASA1
About Gene Symbol	RASA1 (RAS P21 Protein Activator 1) is a GTPase-activating protein involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of RASA1 in A549, HL-60, U87-MG lysates using RASA1 antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using RASA1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.