

S100 alpha 6 Rabbit mAb


ABZ3361

Product Overview


Name	S100 alpha 6 Rabbit mAb
Accession(Primary)	P06703
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 10 kDa; Observed MW: 10 kDa
Recommended Dilution	WB-1:1000-1:5000; ICC/IF-1:50-1:200; FC-1:100-1:1000

Target information (P06703)


Gene Name	S100A6
About Gene Symbol	



Western blot analysis of S100 alpha 6/PRA in 3T3 lysates using S100 alpha6 antibody.



Western blot analysis of S100 alpha 6/PRA in HeLa, A549 lysates using S100 alpha 6/PRA antibody.



Immunocytochemistry analysis of S100 alpha 6/PRA (green) in hela using S100 alpha 6/PRA antibody, and DAPI(blue)