


SARA Rabbit mAb

ABZ3362


Product Overview

Name	SARA Rabbit mAb
Accession(Primary)	O95405
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Rat
Calculated MW	Calculated MW: 156 kDa; Observed MW: 156 kDa
Recommended Dilution	WB-1:2000-1:5000; IHC-P-1:100-1:200

Target information (O95405)

Gene Name ZFYVE9**About Gene Symbol**

Western blot analysis of SARA in HeLa, A549, HL-60 lysates using SARA antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using SARA antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.