

SEC24D Rabbit mAb


ABZ3363

Product Overview


Name	SEC24D Rabbit mAb
Accession(Primary)	O94855
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 113 kDa; Observed MW: 113 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (O94855)

Gene Name	SEC24D
About Gene Symbol	SEC24D (SEC24 Homolog D, COPII Coat Complex Component) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely st



Western blot analysis of SEC24D in Raw264.7, C6 lysates using SEC24D antibody.



Western blot analysis of SEC24D in HeLa, A549, HL-60 lysates using SEC24D antibody.