

Semaphorin 7A Rabbit mAb


ABZ3364

Product Overview


Name	Semaphorin 7A Rabbit mAb
Accession(Primary)	O75326
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 75 kDa; Observed MW: 75 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:50-1:100; FC-1:200-1:500; IP-1:20-1:50

Target information (O75326)

Gene Name	SEMA7A
About Gene Symbol	



Western blot analysis of Semaphorin 7a in A549, HL-60, U251, U87-MG lysates using Semaphorin 7A antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsil using Semaphorin 7a antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.