

SF3B1 Rabbit mAb


ABZ3365

Product Overview


Name	SF3B1 Rabbit mAb
Accession(Primary)	O75533
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 146 kDa; Observed MW: 155 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200; ICC/IF-1:200-1:500; FC-1:200-1:2000

Target information (O75533)


Gene Name	SF3B1
About Gene Symbol	




Immunocytochemistry analysis of SF3B1 (green)
in 293T using SF3B1 antibody, and DAPI (blue).



Western blot analysis of SF3B1 in 3T3 lysates
using SF3B1 antibody.



Western blot analysis of SF3B1 in HeLa, CHO-K1, C6, Ramos lysates using SF3B1 antibody



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using SF3B1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.