

SV40 T Antigen Rabbit mAb


ABZ3367

Product Overview


Name	SV40 T Antigen Rabbit mAb
Accession(Primary)	P03070
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Virus
Calculated MW	Calculated MW: 82 kDa; Observed MW: 90 kDa
Recommended Dilution	WB-1:500-1:1000; IF-1:20-1:50; FC-1:200-1:500; IP-1:20-1:50

Target information (P03070)

Gene Name**About Gene Symbol**



Western blot analysis of SV40 T antigen in A549, HL-60, U251 lysates using SV40 T Antigen antibody.



Immunocytochemistry analysis of SV40 T antigen (green) in HeLa using SV40 T antigen antibody, and DAPI (blue)