

TATA binding protein TBP Rabbit mAb


ABZ3368

Product Overview


Name	TATA binding protein TBP Rabbit mAb
Accession(Primary)	P20226
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 38 kDa; Observed MW: 38 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:500-1:2000; ICC/IF-1:200-1:500; FC-1:200-1:500; IP-1:20-1:50

Target information (P20226)


Gene Name	TBP
About Gene Symbol	




Western blot analysis of TATA binding protein TBP in 3T3 lysates using TATA Box Binding Protein antibody.



Western blot analysis of TATA binding protein TBP in HeLa, A549, HL-60, U2OS lysates using TATA binding protein TBP antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsil using TATA binding protein TBP antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunocytochemistry analysis of TATA binding protein TBP (green) in HeLa using TATA binding protein TBP antibody, and DAPI (blue).