

TCP 1 alpha Rabbit mAb



ABZ3369

Product Overview

Name	TCP 1 alpha Rabbit mAb
Accession(Primary)	P17987
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 60 kDa; Observed MW: 60 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200; IP-1:10-1:100; ICC/IF-1:50-1:100


Target information (P17987)

Gene Name	TCP1
About Gene Symbol	TCP1 (T-Complex 1) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Immunocytochemistry analysis of TCP 1 alpha (green) in HeLa using TCP 1 alpha antibody, and DAPI (blue).

Western blot analysis of TCP1 alpha/CCTA in HeLa, CHO-K1, Ramos lysates using TCP1 alpha/CCTA antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer tissue using TCP1 alpha/CCTA antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.