

TFEB Rabbit mAb


ABZ3370

Product Overview


Name	TFEB Rabbit mAb
Accession(Primary)	P19484
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 53 kDa; Observed MW: 65-70 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:200-1:500; ICC/IF-1:20-1:50; FC-1:200-1:500; IP-1:20-1:50

Target information (P19484)


Gene Name	TFEB
About Gene Symbol	



Immunocytochemistry analysis of Transcription Factor EB (green) in Hela using Transcription Factor EB antibody, and DAPI (blue).



Western blot analysis of TFEB in A549, HL-60 lysates using TFEB antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer tissue using TFEB antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.