

Tissue Factor Pathway Inhibitor Rabbit mAb

ABZ3371


Product Overview

Name	Tissue Factor Pathway Inhibitor Rabbit mAb
Accession(Primary)	P10646
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 35 kDa; Observed MW: 40-55 kDa
Recommended Dilution	WB-1:500-1:1000; IP-1:20

Target information (P10646)

Gene Name

TFPI

About Gene Symbol

Western blot analysis of TFPI in A549, HL-60, U2OS lysates using TFPI antibody.

