

Thrombospondin 1 Rabbit mAb


ABZ3373

Product Overview


Name	Thrombospondin 1 Rabbit mAb
Accession(Primary)	P07996
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 129 kDa; Observed MW: 170 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:2000; ICC/IF-1:50-1:100; FC-1:50-1:500; IP-1:20-1:50

Target information (P07996)

Gene Name	THBS1
About Gene Symbol	



Western blot analysis of Thrombospondin 1 in MCF-7, 293T, C2C12 lysates using Thrombospondin 1 antibody.



Immunohistochemistry analysis of paraffin-embedded Human Cholangiocarcinoma using Thrombospondin 1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.