

TPMT Rabbit mAb


ABZ3375

Product Overview


Name	TPMT Rabbit mAb
Accession(Primary)	P51580
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 28 kDa; Observed MW: 28 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200; FC-1:100-1:500

Target information (P51580)

Gene Name	TPMT
About Gene Symbol	



Immunohistochemistry analysis of paraffin-embedded Human Cholangiocarcinoma using TPMT antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Western blot analysis of TPMT in HeLa, A549, HL-60 lysates using TPMT antibody.