

TRAPPC2 Rabbit mAb

ABZ3377


Product Overview

Name	TRAPPC2 Rabbit mAb
Accession(Primary)	P0DI81
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 16 kDa; Observed MW: 16 kDa
Recommended Dilution	WB-1:1000-1:5000

Target information (P0DI81)

Gene Name TRAPPC2

About Gene Symbol



Western blot analysis of TRAPPC2 in HeLa, HL-60 lysates using TRAPPC2 antibody.

