

TRIM24 Rabbit mAb


ABZ3378

Product Overview

Name	TRIM24 Rabbit mAb
Accession(Primary)	O15164
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 117 kDa; Observed MW: 117 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:100; IF-1:50-1:200

Target information (O15164)

Gene Name	TRIM24
About Gene Symbol	TRIM24 (Tripartite Motif Containing 24) is a ubiquitin E3 ligase involved in ubiquitin-proteasome pathway. It is an important research target in gene regulation, development, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of TRIM24 in HeLa, A549
lysates using TRIM24 antibody.