

Ubiquitin Conjugating Enzyme E2 C Rabbit mAb


ABZ3379

Product Overview


Name	Ubiquitin Conjugating Enzyme E2 C Rabbit mAb
Accession(Primary)	O00762
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 20 kDa; Observed MW: 20 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:200-1:500; ICC/IF-1:20-1:50; IP-1:20-1:50

Target information (O00762)


Gene Name	UBE2C
About Gene Symbol	




Immunocytochemistry analysis of Ubiquitin Conjugating Enzyme E2 C (green) in HeLa using Ubiquitin Conjugating Enzyme E2 C antibody, and DAPI (blue).



Western blot analysis of UBE2C in PC-12, Raw264.7, 3T3 lysates using UBE2C antibody.



Western blot analysis of UBE2C in A549, HL-60, U2OS lysates using UBE2C antibody



Immunohistochemistry analysis of paraffin-embedded Human tonsil using UBE2C antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.