

Ubiquitin K63 Rabbit mAb


ABZ3381

Product Overview


Name	Ubiquitin K63 Rabbit mAb
Accession(Primary)	P0CG47
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 26 kDa; Observed MW: 8 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:100; IF-1:50-1:200; ICC/IF-1:50-1:100

Target information (P0CG47)


Gene Name	UBB
About Gene Symbol	



Western blot analysis of Ubiquitin (linkagespecific K63) in HeLa, A549 lysates using Ubiquitin K63 antibody.



Immunohistochemistry analysis of paraffin-embedded Human colon cancer using Ubiquitin (linkagespecific K63) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunocytochemistry analysis of Ubiquitin
(linkagespecific K63) (green) in hela using Ubiquitin
(linkagespecific K63) antibody, and DAPI(blue)