

## XRCC1 Rabbit mAb


ABZ3383

### Product Overview


|                             |                                                            |
|-----------------------------|------------------------------------------------------------|
| <b>Name</b>                 | XRCC1 Rabbit mAb                                           |
| <b>Accession(Primary)</b>   | P18887                                                     |
| <b>Host Species</b>         | Rabbit                                                     |
| <b>Clonality</b>            | Monoclonal                                                 |
| <b>Isotype</b>              | IgG                                                        |
| <b>Species Reactivity</b>   | Human, Mouse, Rat                                          |
| <b>Calculated MW</b>        | Calculated MW: 69 kDa; Observed MW: 82 kDa                 |
| <b>Recommended Dilution</b> | WB-1:1000-1:5000; IHC-P-1:200-1:500;<br>ICC/IF-1:100-1:200 |

### Target information (P18887)


|                          |       |
|--------------------------|-------|
| <b>Gene Name</b>         | XRCC1 |
| <b>About Gene Symbol</b> |       |



Western blot analysis of XRCC1 in PC-12 lysates using XRCC1 antibody.



Western blot analysis of XRCC1 in A549, MCF-7, 293T lysates using XRCC1 antibody.



Immunocytochemistry analysis of XRCC1 (green)  
in hela using XRCC1 antibody, and DAPI(blue)