

ZBTB7A Rabbit mAb


ABZ3385

Product Overview


Name	ZBTB7A Rabbit mAb
Accession(Primary)	O95365
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 61 kDa; Observed MW: 61 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200; ICC/IF-1:200-1:500; IP-1:10-1:100

Target information (O95365)


Gene Name	ZBTB7A
About Gene Symbol	



Western blot analysis of ZBTB7A in HeLa, Ramos lysates using ZBTB7A antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsil using ZBTB7A antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunocytochemistry analysis of ZBTB7A (green)
in MCF-7 using ZBTB7A antibody, and DAPI (blue)