


DAPK3 Rabbit mAb

ABZ3386

Product Overview

Name	DAPK3 Rabbit mAb
Accession(Primary)	O43293
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 53 kDa; Observed MW: 53 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (O43293)

Gene Name DAPK3**About Gene Symbol**

Western blot analysis of ZIP Kinase in HeLa, A549, MCF-7, U251 lysates using DAP kinase 3 antibody.

