

PIAS1/2 Rabbit mAb


ABZ3388

Product Overview


Name	PIAS1/2 Rabbit mAb
Accession(Primary)	O75925
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 72 kDa; Observed MW: 76 kDa
Recommended Dilution	WB-1:500-1:1000; IF-1:100-1:200; ICC/IF-1:200-1:500; IHC-1:50-1:100

Target information (O75925)

Gene Name	PIAS1/PIAS2/PIAS3
About Gene Symbol	



Western blot analysis of PIAS1/PIAS2 in mouse testis, rat testis lysates using PIAS1/2 antibody.



Western blot analysis of PIAS1+PIAS2 in HeLa, A549, HL-60, U2OS lysates using PIAS1+PIAS2 antibody.