

MEF2A/MEF2C Rabbit mAb


ABZ3390

Product Overview


Name	MEF2A/MEF2C Rabbit mAb
Accession(Primary)	Q06413
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 51 kDa; Observed MW: 51 kDa
Recommended Dilution	WB-1:500-1:1000; IHC:1:100-1:200; ICC/IF-1:100-1:200; FC-1:100-1:200

Target information (Q06413)

Gene Name	MEF2C/MEF2A
About Gene Symbol	



Western blot analysis of MEF2A+MEF2C in C6, mouse brain, rat brain lysates using MEF2A/MEF2C antibody.



Western blot analysis of MEF2A+MEF2C in HL-60, U2OS, U87-MG lysates using MEF2A+MEF2C antibody.