

Human IgG4 Rabbit mAb

ABZ3391


Product Overview

Name	Human IgG4 Rabbit mAb
Accession(Primary)	P01861
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 36 kDa; Observed MW: 36 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:500-1:5000


Target information (P01861)

Gene Name IGHG4

About Gene Symbol



Western blot analysis of IgG4 in HL-60, U2OS lysates using Human IgG4 antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsil using IgG4 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.