

Human kappa Light Chain Rabbit mAb


ABZ3394

Product Overview


Name	Human kappa Light Chain Rabbit mAb
Accession(Primary)	P01834
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 12 kDa; Observed MW: 26 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:500; ICC/IF-1:200-1:500

Target information (P01834)


Gene Name	IGKC
About Gene Symbol	



Immunocytochemistry analysis of Human kappa Light Chain (green) in HeLa using Human kappa Light Chain antibody, and DAPI (blue).



Western blot analysis of Human Kappa Light Chain in A549, HL-60, U2OS lysates using Human Kappa Light Chain antibody. Observed band size: 26 kDa.



Immunohistochemistry analysis of paraffin-embedded Human tonsil tissue using Human Kappa Light Chain antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.