

PRPF31 Rabbit mAb


ABZ3396

Product Overview


Name	PRPF31 Rabbit mAb
Accession(Primary)	Q8WWY3
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 55 kDa; Observed MW: 55 kDa
Recommended Dilution	WB-1:1000-1:5000; ICC/IF-1:100-1:200; FC-1:100-1:200

Target information (Q8WWY3)


Gene Name	PRPF31
About Gene Symbol	



Western blot analysis of PRPF31 in PC-12 lysates using PRPF31 antibody.



Western blot analysis of PRPF31 in C2C12, PC-3 lysates using PRPF31 antibody.



Immunocytochemistry analysis of PRPF31 (green)
in HeLa using PRPF31 antibody, and DAPI (blue)