

CD147 Rabbit mAb


ABZ3397

Product Overview


Name	CD147 Rabbit mAb
Accession(Primary)	P18572
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Mouse
Calculated MW	Calculated MW: 42 kDa; Observed MW: 38-58 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:1000-1:5000; ICC/IF-1:100-1:200; FC-1:100-1:200; IP-1:20-1:50

Target information (P18572)

Gene Name	Bsg
About Gene Symbol	



Immunocytochemistry analysis of CD147 (green) in 3T3 using CD147 antibody, and DAPI (blue).



Western blot analysis of CD147 in rat Brain, C6, 3T3 lysates using CD147 antibody.