


Collagen III alpha 1 Rabbit mAb

ABZ3398


Product Overview

Name	Collagen III alpha 1 Rabbit mAb
Accession(Primary)	P02461
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 139 kDa; Observed MW: 150 kDa
Recommended Dilution	WB-1:1000-1:4000; IF: 1:50; IP-1:20

Target information (P02461)

Gene Name COL3A1**About Gene Symbol**

Western blot analysis of Collagen III in K562, C6, 3T3 lysates using Collagen III alpha 1 antibody.



Immunocytochemistry analysis of Collagen
III(green) in Hela using Collagen III antibody, and
DAPI(blue)