


## S6K1 Rabbit mAb


ABZ3399

### Product Overview

|                             |                                            |
|-----------------------------|--------------------------------------------|
| <b>Name</b>                 | S6K1 Rabbit mAb                            |
| <b>Accession(Primary)</b>   | P23443                                     |
| <b>Host Species</b>         | Rabbit                                     |
| <b>Clonality</b>            | Monoclonal                                 |
| <b>Isotype</b>              | IgG                                        |
| <b>Species Reactivity</b>   | Human, Mouse, Rat                          |
| <b>Calculated MW</b>        | Calculated MW: 59 kDa; Observed MW: 70 kDa |
| <b>Recommended Dilution</b> | WB-1:1000-1:5000; IHC-P-1:20-1:50          |

### Target information (P23443)

**Gene Name** RPS6KB1**About Gene Symbol**



Western blot analysis of S6K1 in 293 lysates using  
S6K1 antibody.