

Vitamin D Binding protein Rabbit mAb


ABZ3400

Product Overview


Name	Vitamin D Binding protein Rabbit mAb
Accession(Primary)	P02774
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 53 kDa; Observed MW: 53 kDa
Recommended Dilution	WB-1:1000-1:5000; FC-1:20-1:50; ICC/IF-1:50-1:100

Target information (P02774)

Gene Name	GC
About Gene Symbol	



Immunocytochemistry analysis of Vitamin D Binding Protein (green) in HEPG2 using Vitamin D Binding Protein antibody, and DAPI (blue).



Western blot analysis of Vitamin D Binding protein in rat Brain, 3T3, Hela lysates using Vitamin D Binding protein antibody.