

Klotho Rabbit mAb

ABZ3403


Product Overview

Name	Klotho Rabbit mAb
Accession(Primary)	Q9UEF7
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 116 kDa; Observed MW: 62 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:200-1:500


Target information (Q9UEF7)

Gene Name KL

About Gene Symbol



Western blot analysis of Klotho in Human fetal kidney lysates using Klotho antibody.



Immunohistochemistry analysis of paraffin-embedded Human fetal kidney using klotho antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.