

## Bcl 2 Rabbit mAb


ABZ3404

### Product Overview


|                             |                                                     |
|-----------------------------|-----------------------------------------------------|
| <b>Name</b>                 | Bcl 2 Rabbit mAb                                    |
| <b>Accession(Primary)</b>   | P10415                                              |
| <b>Host Species</b>         | Rabbit                                              |
| <b>Clonality</b>            | Monoclonal                                          |
| <b>Isotype</b>              | IgG                                                 |
| <b>Species Reactivity</b>   | Human                                               |
| <b>Calculated MW</b>        | Calculated MW: 26 kDa; Observed MW: 26 kDa          |
| <b>Recommended Dilution</b> | WB-1:500-1:1000; IHC-P-1:200-1:500;<br>IP-1:20-1:50 |

### Target information (P10415)

|                          |                                                                                                                                                                                                                                |
|--------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Gene Name</b>         | BCL2                                                                                                                                                                                                                           |
| <b>About Gene Symbol</b> | BCL2 (BCL2 Apoptosis Regulator) is a protein phosphatase involved in apoptosis. It is an important research target in cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets. |



Western blot analysis of Bcl2 in Jurkat, HeLa lysates using Bcl2 antibody.



Immunohistochemistry analysis of paraffin-embedded Human breast cancer using Bcl2 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.