


N Cadherin Rabbit mAb

ABZ3408


Product Overview

Name	N Cadherin Rabbit mAb
Accession(Primary)	P19022
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 100 kDa; Observed MW: 140 kDa
Recommended Dilution	WB-1:500-1:1000; IP-1:20

Target information (P19022)

Gene Name CDH2**About Gene Symbol**

Western blot analysis of N Cadherin in mouse brain, rat brain lysates using N Cadherin antibody.



Western blot analysis of N Cadherin in HepG2
lysates using N Cadherin antibody.