

Cathepsin B Rabbit mAb


ABZ3410

Product Overview


Name	Cathepsin B Rabbit mAb
Accession(Primary)	P10605
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Mouse, Rat
Calculated MW	Calculated MW: 37 kDa; Observed MW: 44,27,24 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:200-1:1000; ICC/IF-1:100-1:200; FC-1:200-1:500

Target information (P10605)

Gene Name	Ctsb
About Gene Symbol	



Western blot analysis of Cathepsin B in rat Brain, C6, 3T3 lysates using Cathepsin B antibody.



Immunocytochemistry analysis of Cathepsin B (green) in Raw264.7 using Cathepsin B antibody, and DAPI (blue)