

FACL4 Rabbit mAb


ABZ3411

Product Overview


Name	FACL4 Rabbit mAb
Accession(Primary)	O60488
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 79 kDa; Observed MW: 79 kDa
Recommended Dilution	WB-1:2000-1:10000; IHC-1:200-1:1000; ICC/IF-1:100-1:200; FC-1:10-1:100; IP-1:10-1:100

Target information (O60488)

Gene Name	ACSL4
About Gene Symbol	ACSL4 (Acyl-CoA Synthetase Long Chain Family Member 4) is a biosynthetic enzyme involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of FACL4 in K562, C6, 3T3, HeLa lysates using FACL4 antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using FACL4 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.