

Heme Oxygenase 1 Rabbit mAb


ABZ3412

Product Overview


Name	Heme Oxygenase 1 Rabbit mAb
Accession(Primary)	P09601
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 33 kDa; Observed MW: 33 kDa
Recommended Dilution	WB-1:2000-1:10000; IHC-1:100-1:400; IF-1:50-1:200; FC-1:20; IP-1:20

Target information (P09601)


Gene Name	HMOX1
About Gene Symbol	




Western blot analysis of Heme Oxygenase 1 in rat spleen lysates using Heme Oxygenase 1 antibody.



Western blot analysis of Heme Oxygenase 1 in mouse spleen lysates using Heme Oxygenase 1 antibody.



Western blot analysis of Heme Oxygenase 1 in U251MG, K562, A549 lysates using Heme Oxygenase 1 antibody



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using Heme Oxygenase 1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.