

Iba1 Rabbit mAb


ABZ3413

Product Overview


Name	Iba1 Rabbit mAb
Accession(Primary)	O70200
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 17 kDa; Observed MW: 17 kDa
Recommended Dilution	WB-1:500-1:1000; IHC:1:200-1:1000; ICC/IF-1:100-1:1000; FC-1:50-1:100

Target information (O70200)

Gene Name	Aif1
About Gene Symbol	



Western blot analysis of Iba1 in C6, 3T3, HeLa lysates using Iba1 antibody.



Immunohistochemistry analysis of paraffin-embedded Human colon cancer using Iba1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.