


Liver Arginase Rabbit mAb

ABZ3415


Product Overview

Name	Liver Arginase Rabbit mAb
Accession(Primary)	P07824
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Mouse, Rat
Calculated MW	Calculated MW: 35 kDa; Observed MW: 35 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-:1:200-1:2000

Target information (P07824)

Gene Name Arg1**About Gene Symbol**

Western blot analysis of Liver Arginase in rat liver, PC-12 lysates using Liver Arginase antibody.



Western blot analysis of Liver Arginase in mouse liver lysates using Liver Arginase antibody.